

UNIVERSITY OF CAMBRIDGE INTERNATIONAL EXAMINATIONS General Certificate of Education Advanced Subsidiary Level and Advanced Level

	CANDIDATE NAME		
	CENTRE NUMBER	CANDIDATE NUMBER	
* 1 8	CHEMISTRY		9701/34
6 u	Advanced Prac	tical Skills O	ctober/November 2010
ω ω O 7 O 20 5	Candidates ans Additional Mate	wer on the Question Paper. rials: As listed in the Instructions to Supervisors	
*	READ THESE I	NSTRUCTIONS FIRST	
	Give details of t Write in dark bl You may use a Do not use stap	tre number, candidate number and name on all the work you hand in. he practical session and laboratory where appropriate, in the boxes p ue or black pen. soft pencil for any diagrams, graphs or rough working. bles, paper clips, highlighters, glue or correction fluid. E IN ANY BARCODES.	rovided.
		arks if you do not show your working or if you do not use	
	appropriate unit Use of a Data E	Booklet is unnecessary.	Session
Qualitative Analysis Notes are printed on pages 13 and 14.			
		e examination, fasten all your work securely together. marks is given in brackets [] at the end of each question or	Laboratory

For Examiner's Use	
1	
2	
3	
Total	

This document consists of **12** printed pages and **4** blank pages.



There are three questions on this paper. Question 2 should not be the last question attempted.

FB 1 is an aqueous solution containing 21.50 gdm⁻³ of a mixture of iron(II) sulfate, FeSO₄ and iron(III) sulfate, Fe₂(SO₄)₃.
FB 2 is an aqueous solution containing 2.00 gdm⁻³ potassium manganate(VII), KMnO₄.

In the presence of acid, the iron(II) sulfate is oxidised by potassium manganate(VII).

 $2KMnO_4(aq) + 8H_2SO_4(aq) + 10FeSO_4(aq) \rightarrow 5Fe_2(SO_4)_3(aq) + 2MnSO_4(aq) + K_2SO_4(aq) + 8H_2O(l)$

(a) Method

- Fill a burette with **FB 2**.
- Pipette 25.0 cm³ of **FB 1** into the conical flask.
- Use a 25 cm³ measuring cylinder to add 10 cm³ of dilute sulfuric acid to the flask.
- Place the flask on a white tile.
- Carefully titrate with **FB 2** until the first permanent pink colour is obtained.

You should perform a **rough titration**.

In the space below record your burette readings for this rough titration.

The rough titre is cm³.

- Carry out as many accurate titrations as you think necessary to obtain consistent results.
- Record in a suitable form below all of your burette readings and the volume of **FB 2** added in each accurate titration.
- Make certain any recorded results show the precision of your practical work.

For

Examiner's Use

[7]

(b) From your accurate titration results obtain a suitable value to be used in your calculation. Show clearly how you have obtained this value.

3

For Examiner's Use

25.0 cm³ of **FB 1** required cm³ of **FB 2**. [1]

Calculations

Show your working and appropriate significant figures in the final answer to **each** step of your calculations.

(c) (i) Calculate the concentration, in moldm⁻³, of the potassium manganate(VII) in FB 2.
FB 2 contains 2.00 g dm⁻³ KMnO₄.
[A_r: O, 16.0; K, 39.1; Mn, 54.9]

The concentration of potassium manganate(VII) in **FB 2** is moldm⁻³.

(ii) Calculate how many moles of $KMnO_4$ were present in the volume calculated in (b).

..... mol of $KMnO_4$.

(iii) Calculate how many moles of iron(II) sulfate, FeSO₄, reacted with the potassium manganate(VII) in (ii).

 $2\mathsf{KMnO}_4(\mathsf{aq}) + 8\mathsf{H}_2\mathsf{SO}_4(\mathsf{aq}) + 10\mathsf{FeSO}_4(\mathsf{aq}) \longrightarrow 5\mathsf{Fe}_2(\mathsf{SO}_4)_3(\mathsf{aq}) + 2\mathsf{MnSO}_4(\mathsf{aq}) + \mathsf{K}_2\mathsf{SO}_4(\mathsf{aq}) + 8\mathsf{H}_2\mathsf{O}(\mathsf{I})$

..... mol of $FeSO_4$ reacted with the potassium manganate(VII).

(iv)	Calculate the concentration, in moldm ^{-3} of FeSO ₄ in FB 1 .	For Examiner's Use
	The concentration of $FeSO_4$ in FB 1 is mol dm ⁻³ .	
(v)	Calculate the concentration, in g dm ⁻³ , of FeSO ₄ in FB 1 . [A_r : O, 16.0; S, 32.1; Fe, 55.8]	
	FB 1 contains gdm ⁻³ of FeSO₄.	
(vi)	FB 1 is an aqueous solution containing 21.50 g dm ⁻³ of FeSO ₄ and Fe ₂ (SO ₄) ₃ . Calculate the percentage, by mass, of FeSO ₄ in this mixture.	
	The mixture contains % FeSO ₄ . [5]	
	[Total: 13]	

4

2 FB 3 is a mixture containing anhydrous sodium carbonate, Na₂CO₃, and sodium hydrogencarbonate, NaHCO₃.

When heated, sodium hydrogencarbonate decomposes.

 $2NaHCO_3(s) \rightarrow Na_2CO_3(s) + CO_2(g) + H_2O(g)$

Anhydrous sodium carbonate does not decompose when heated.

You are to determine if sodium hydrogencarbonate is the major component, by mass, of the mixture in **FB 3**.

(a) Method – Read through the instructions before starting any practical work.

- Weigh and record the mass of an empty boiling-tube.
- Tip the contents of the tube labelled **FB 3** into the weighed boiling-tube. Reweigh and record the mass of the boiling-tube and **FB 3**.
- Gently heat the **FB 3** in the boiling-tube for 2 minutes then heat strongly for a further 2 minutes.
 - Take care not to lose any solid from the tube during heating.
- Warm the upper parts of the boiling-tube to evaporate any water that may have condensed while heating the solid.
- Place the hot tube on a heat-proof mat and leave to cool.
- You are advised to continue with part (d) of this question or to start another question while the tube cools.
- When cool, reweigh the boiling-tube and the residual sodium carbonate.
- Reheat, cool and reweigh the tube until you are satisfied decomposition is complete.

Results

In an appropriate form, in the space below, record all of your balance readings, the mass of **FB 3** heated, the mass of residual sodium carbonate and the mass loss on heating.





..... [2]

(c)				For
	Suggest to on heating		show that sodium carbonate does not decompose	Examiner's Use
	(i)			
	(ii)			
			[2]	
(d)	A student	is asked to weigh, with max	imum precision, a solid.	
	The three	balances available are:		
	balan	ice A, reading to 1 decimal p ice B, reading to 2 decimal p ice C, reading to 3 decimal p	blaces,	
		llest division on a burette is imum error in a single buret		
	Balance	readings can be treated in t	he same way.	
	Complete	the following table.		
ł	balance maximum error for a single balance reading / g maximum % error when weighing:			
	А	±	9.0g of solid =	
	В	±	4.00 g of solid =	
	С	±	0.500 g of solid =	

[2]

[Total: 12]

3 FB 4, **FB 5**, **FB 6** and **FB 7** are aqueous solutions each containing one of the ions Al^{3+} , NH_4^{++} , Mg^{2+} , Mn^{2+} .

You will carry out the following tests on each of the solutions.

At each stage of any test you are to record details of the following.

- colour changes seen
- the formation of any precipitate
- the solubility of such precipitates in an excess of the reagent added

Where gases are released they should be identified by a test, **described in the appropriate place in your observations**.

You should indicate clearly at what stage in a test a change occurs. Marks are **not** given for chemical equations. **No additional tests for ions present should be attempted.** If any solution is warmed, a boiling-tube MUST be used.

(a) Carry out the following tests. Record your observations in the spaces provided in the table.

toot		observations				
	test	FB 4	FB 5	FB 6	FB 7	
(i)	To 1 cm depth of solution in a test- tube add 1 cm depth of aqueous sodium hydroxide. Swirl the tube, then					
	add a further 2 cm depth of aqueous sodium hydroxide.					
	In tests (ii) and (iii) µ	out a cross in a	ny boxes where	the test is not ca	arried out.	
(ii)	If a precipitate remains at the end of test (i) leave the test-tube and contents to stand for a few minutes.					
(iii)	If no precipitate formed at all in test (i) tip the contents of the tube into a boiling-tube and warm gently. Care: heated solutions containing sodium hydroxide are liable to be ejected from the tube.					

© UCLES 2010

9701/34/O/N/10

For Examiner's Use

test		observations				
		FB 4	FB 5	FB 6	FB 7	Use
(iv)	To 1 cm depth of solution in a test- tube add 1 cm depth of aqueous ammonia. Swirl the tube, then					
	add a further 2 cm depth of aqueous ammonia.					V

11

[6]

I

Ш

Ш

IV

(b) Use the Qualitative Analysis Notes on page 13 to identify the cation present in each of the solutions.

Complete the table below to identify each ion and to give supporting evidence from your observations.

solution	cation	supporting evidence
FB 4		
FB 5		
FB 6		
FB 7		

[4]

Rinse and re-use test-tubes where possible.

(c) Carry out the following tests on the solution you have identified as containing Al^{3+} ions and record your observations in the spaces provided.

For Examiner's Use

		observation
(i)	Add aqueous sodium iodide	
(ii)	Add dilute sulfuric acid	

Explain how your results confirm the presence of Al^{3+} and eliminate any other ion.

(d) What other cation listed in the Qualitative Analysis Notes on page 13 would give similar results to Al³⁺ in (a)?

.....[1]

(e) Carry out the following tests and make careful observations of all that happens in each experiment. Complete the table.

	test	observations
(i)	To 1 cm depth of aqueous silver nitrate in a test-tube add 1 cm depth of aqueous sodium chloride. Keep the tube for comparison with the observations in test (ii) .	
(ii)	Repeat test (i) . To 1 cm depth of aqueous silver nitrate in a test-tube add 1 cm depth of aqueous sodium chloride, then	Do not repeat your observations from test (i)
	add 1 cm depth of aqueous sodium iodide and shake the tube.	

(f) Suggest an explanation for your observations when aqueous sodium iodide is added in test (e)(ii).

[1]

[2]

Qualitative Analysis Notes

Key: [ppt. = precipitate]

1 Reactions of aqueous cations

	reaction with				
ion	NaOH(aq)	NH ₃ (aq)			
aluminium, A <i>l</i> ³⁺ (aq)	white ppt. soluble in excess	white ppt. insoluble in excess			
ammonium, NH ₄ +(aq)	no ppt. ammonia produced on heating				
barium, Ba ²⁺ (aq)	no ppt. (if reagents are pure)	no ppt.			
calcium, Ca ²⁺ (aq)	white ppt. with high [Ca ²⁺ (aq)]	no ppt.			
chromium(III), Cr ³⁺ (aq)	grey-green ppt. soluble in excess giving dark green solution	grey-green ppt. insoluble in excess			
copper(II), Cu ²⁺ (aq)	pale blue ppt. insoluble in excess	blue ppt. soluble in excess giving dark blue solution			
iron(II), Fe ²⁺ (aq)	green ppt. turning brown on contact with air insoluble in excess	green ppt. turning brown on contact with air insoluble in excess			
iron(III), Fe ³⁺ (aq)	red-brown ppt. insoluble in excess	red-brown ppt. insoluble in excess			
lead(II), Pb ²⁺ (aq)	white ppt. soluble in excess	white ppt. insoluble in excess			
magnesium, Mg ²⁺ (aq)	white ppt. insoluble in excess	white ppt. insoluble in excess			
manganese(II), Mn ²⁺ (aq)	off-white ppt. rapidly turning brown on contact with air insoluble in excess	off-white ppt. rapidly turning brown on contact with air insoluble in excess			
zinc, Zn ²⁺ (aq)	white ppt. soluble in excess	white ppt. soluble in excess			

[Lead(II) ions can be distinguished from aluminium ions by the insolubility of lead(II) chloride.]

2 Reactions of anions

ion	reaction
carbonate,	CO ₂ liberated by dilute acids
CO3 ²⁻	
chromate(VI), CrO ₄ ²⁻ (aq)	yellow solution turns orange with H ⁺ (aq); gives yellow ppt. with Ba ²⁺ (aq); gives bright yellow ppt. with Pb ²⁺ (aq)
chloride, C <i>l</i> ⁻ (aq)	gives white ppt. with Ag ⁺ (aq) (soluble in NH ₃ (aq)); gives white ppt. with Pb ²⁺ (aq)
bromide, Br ⁻ (aq)	gives white ppt. with Ag ⁺ (aq) (partially soluble in NH ₃ (aq)); gives white ppt. with Pb ²⁺ (aq)
iodide, I⁻ (aq)	gives yellow ppt. with Ag ⁺ (aq) (insoluble in NH ₃ (aq)); gives yellow ppt. with Pb ²⁺ (aq)
nitrate, NO ₃ ⁻ (aq)	NH_3 liberated on heating with OH ⁻ (aq) and A <i>l</i> foil
nitrite, NO ₂ ⁻ (aq)	NH ₃ liberated on heating with OH [−] (aq) and A <i>l</i> foil; NO liberated by dilute acids (colourless NO → (pale) brown NO ₂ in air)
sulfate, $SO_4^{2^-}$ (aq)	gives white ppt. with Ba ²⁺ (aq) or with Pb ²⁺ (aq) (insoluble in excess dilute strong acid)
sulfite, SO ₃ ²⁻ (aq)	SO ₂ liberated with dilute acids; gives white ppt. with Ba ²⁺ (aq) (soluble in excess dilute strong acid)

3 Tests for gases

gas	test and test result
ammonia, NH ₃ turns damp red litmus paper blue	
carbon dioxide, CO ₂	gives a white ppt. with limewater (ppt. dissolves with excess CO ₂)
chlorine, Cl ₂	bleaches damp litmus paper
hydrogen, H ₂	"pops" with a lighted splint
oxygen, O ₂	relights a glowing splint
sulfur dioxide, SO ₂	turns acidified aqueous potassium dichromate(VI) from orange to green

15

16

Permission to reproduce items where third-party owned material protected by copyright is included has been sought and cleared where possible. Every reasonable effort has been made by the publisher (UCLES) to trace copyright holders, but if any items requiring clearance have unwittingly been included, the publisher will be pleased to make amends at the earliest possible opportunity.

University of Cambridge International Examinations is part of the Cambridge Assessment Group. Cambridge Assessment is the brand name of University of Cambridge Local Examinations Syndicate (UCLES), which is itself a department of the University of Cambridge.